




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Alternator Vacuum Pump Cutaway

**Product Code :**  
LBNY-0005-121000408



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Alternator Vacuum Pump Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

**Technical Specification :**

Operated: Manually  
Dimensions: 20 x 20 x 25(h) cm  
Net Weight: 8 kg  
Gross Weight: 15 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**