



Email : sales@laboratoryinstrumentindia.com

Product Name :
Alternator Liquid Cooled Windingless Rotor Compact
Cutaway

Product Code :
LBNY-0005-121000409



Description :

Many parts have been chromium-plated and galvanized for a longer life. Alternator Liquid Cooled Windingless Rotor Compact Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

Technical Specification :

Operated: Manually
Dimensions: 20 x 20 x 25(h) cm
Net Weight: 8 kg
Gross Weight: 15 kg
Equipment Layout: On Base.



Laboratory Instrument India