




Email : sales@laboratoryinstrumentindia.com

Product Name :
Two Stroke Direct Injection Diesel Engine Model Cutaway

Product Code :
LBNY-0005-121000410



Description :

Many parts have been chromium-plated and galvanized for a longer life. Two Stroke Direct Injection Diesel Engine Model Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 40 x 40 x 65(h) cm
Net Weight: 9 kg
Gross Weight: 13 kg
Equipment Layout: On Base.



Laboratory Instrument India