



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Flow Sensor Cutaway

Product Code :
LBNY-0005-121000411



Description :

Many parts have been chromium-plated and galvanized for a longer life. Air Flow Sensor Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

Technical Specification :

Operated: Static
Dimensions: 10 x 10 x 10(h) cm
Net Weight: 1 kg
Gross Weight: 1.5 kg
Equipment Layout: On Base.



Laboratory Instrument India