



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Air Flow Sensor Cutaway

**Product Code :**  
LBNY-0005-121000411



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Air Flow Sensor Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

**Technical Specification :**

Operated: Static  
Dimensions: 10 x 10 x 10(h) cm  
Net Weight: 1 kg  
Gross Weight: 1.5 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**