




Email : sales@laboratoryinstrumentindia.com

Product Name :
Opposed Piston Engine Cutaway

Product Code :
LBNY-0005-121000413



Description :

Many parts have been chromium-plated and galvanized for a longer life. Opposed Piston Engine Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

Air cooling
Ignition with magneto
Gear distribution with camshaft in the crankcase
Single-body carburettor.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 140 x 110 x 150(h) cm
Net Weight: 190 kg
Gross Weight: 280 kg
Equipment Layout: On Stand with Wheels.

Laboratory Instrument India

