



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Soap Stone

Product Code :
LBNY-0004-740010



Description :

Workshop Soap Stone

Technical Specification :

Soapstone

White

Flat or round soap stones with flat aluminum holder

flat and round holder with soapstone chalk, 4pcs flat and 4pcs round soapstone replacement.

#1 grade soapstone, mark white clear and removable markings.

High temperature resistance and waterproof soapstone, The marks are much clearer and won't be burnt off during welding.

A combination of flat and round aluminum holders, apply to all kinds of operation requirements.

Widely used in welding kit, fabricating project and cutting application. Perfect for marking targets and measurements.

Workshop Soap Stone, Workshop Soap Stone Equipments, Workshop Soap Stone Tools, Workshop Soap Stone Tool Kits, Workshop Soap Stone Manufacturers, Workshop Soap Stone Suppliers from India, China, Kenya

Laboratory Instrument India

