



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbo Jet Engine Model Cutaway

Product Code :
LBNY-0005-121000427



www.laboratoryinstrumentindia.com

Description :

Low-pressure and high-pressure compressor, low-pressure and high-pressure turbine are readily recognisable, as also are the combustion chambers with the injection nozzles and starting plugs. Turbo Jet Engine Model Cutaway of a modern two-wave turbine shows in detail the construction and the operating system of such a motor.

Technical Specification :

Operated: Electrically (Motor 230 V)
Dimensions: 70 x 28 x 45(h) cm
Net Weight: 13 kg
Gross Weight: 30 kg
Equipment Layout: On Base.



Laboratory Instrument India