



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Gases

**Product Code :**  
LBNY-0004-740011



**Description :**

Workshop Gases

**Technical Specification :**

Gases

Argoshield (MAGS)

Argon (TAGS).

Oxygen

Acetylene

Argon (Ar) is present in atmospheric air in very small quantity (about 0.9%). It is colorless, odorless gas and highly inert (non-reactive) like Nitrogen.

Argon is mostly used as an inert shielding gas in welding and other high-temperature industrial processes where ordinarily unreactive substances become reactive.

Workshop Gases, Workshop Gases Equipments, Workshop Gases Tools, Workshop Gases Tool Kits, Workshop Gases Manufacturers, Workshop Gases Suppliers from India, China, Kenya



**Laboratory Instrument India**

