



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motor Cultivator Cutaway

Product Code :
LBNY-0005-121000433



Description :

Many parts have been chromium-plated and galvanized for a longer life. Motor Cultivator Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

Petrol engine 4 strokes single cylinder
Dry clutch with manual control
Handlebars adjustable both in height and width
Motor-stop safety device
Gearbox MTC 3+2
Power: 6 Kw
Recoil starter
Reverse command by lever.

Other Specifications:

Equipment Layout: On Stand with Wheels
Operated: Electrically (Motor 230 V)
Dimensions: 70 x 180 x 110(h) cm
Net Weight: 100 kg
Gross Weight: 160 kg.



Laboratory Instrument India