




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Two Stroke Petrol Engine Model Cutaway

**Product Code :**  
LBNY-0005-121000434



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Two Stroke Petrol Engine Model Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

**Technical Specification :**

Electronic ignition  
Box carburettor  
Piston displacement 46 cc  
Air cooling

**Other Specifications:**

Operated: Manually (Through a Crank Handle)  
Dimensions: 30 x 30 x 40(h) cm  
Net Weight: 8 kg  
Gross Weight: 15 kg  
Equipment Layout: On Base.

---

# Laboratory Instrument India

