




Email : sales@laboratoryinstrumentindia.com

Product Name :
Two Stroke Petrol Engine Model Cutaway

Product Code :
LBNY-0005-121000434



Description :

Many parts have been chromium-plated and galvanized for a longer life. Two Stroke Petrol Engine Model Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

Electronic ignition
Box carburettor
Piston displacement 46 cc
Air cooling

Other Specifications:

Operated: Manually (Through a Crank Handle)
Dimensions: 30 x 30 x 40(h) cm
Net Weight: 8 kg
Gross Weight: 15 kg
Equipment Layout: On Base.

Laboratory Instrument India

