




Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbo Jet Engine Model With High Pressure Compressor
Cutaway

Product Code :
LBNY-0005-121000443



Description :

High-pressure compressor, high pressure turbine are readily recognisable, as well as the combustion chambers with the injection nozzles and starting plugs. Turbo Jet Engine Model With High-Pressure Compressor Cutaway of a modern two-wave turbine shows in detail the construction and the operating system of such a motor.

Technical Specification :

Operated: Electrically (Motor 230 V)
Dimensions: 60 x 40 x 40(h) cm
Net Weight: 10 kg
Gross Weight: 20 kg
Equipment Layout: On Base.



Laboratory Instrument India