




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Turbo Jet Engine Model With High Pressure Compressor  
Cutaway

**Product Code :**  
LBNY-0005-121000443



**Description :**

High-pressure compressor, high pressure turbine are readily recognisable, as well as the combustion chambers with the injection nozzles and starting plugs. Turbo Jet Engine Model With High-Pressure Compressor Cutaway of a modern two-wave turbine shows in detail the construction and the operating system of such a motor.

**Technical Specification :**

Operated: Electrically (Motor 230 V)  
Dimensions: 60 x 40 x 40(h) cm  
Net Weight: 10 kg  
Gross Weight: 20 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**