




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Two Stroke Marine Outboard Engine Cutaway

**Product Code :**  
LBNY-0005-121000449



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Two Stroke Marine Outboard Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

**Technical Specification :**

Water cooling with centrifugal pump  
Mechanical type converter  
2/3 cylinders (2-stroke)  
2-stroke engine

**Other Specifications:**

Operated: Electrically  
Dimensions: 80 x 80 x 150(h) cm  
Net Weight: 60 kg  
Gross Weight: 100 kg  
Equipment Layout: On Stand with Wheels.

---

# Laboratory Instrument India

