




Email : sales@laboratoryinstrumentindia.com

Product Name :
Two Stroke Marine Outboard Engine Cutaway

Product Code :
LBNY-0005-121000449



Description :

Many parts have been chromium-plated and galvanized for a longer life. Two Stroke Marine Outboard Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

Water cooling with centrifugal pump
Mechanical type converter
2/3 cylinders (2-stroke)
2-stroke engine

Other Specifications:

Operated: Electrically
Dimensions: 80 x 80 x 150(h) cm
Net Weight: 60 kg
Gross Weight: 100 kg
Equipment Layout: On Stand with Wheels.

Laboratory Instrument India

