




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Workshop Welding Wire

**Product Code :**  
LBNY-0004-740013



**Description :**

Workshop Welding Wire

**Technical Specification :**

Welding wire

TIG Welding Filler Wire (5) A premium, multiple deoxidized welding filler wire. Produces high-quality welds in most grades of carbon steel. Welding wire is a triple deoxidized wire which provides defect free weld deposits when properly used on most carbon steels. It is used especially for pipe welding.

Welding wire 0.8

TAGS filler wire 1.6 and 2.4 Diameter

Workshop Welding Wire, Workshop Welding Wire Equipments, Workshop Welding Wire Tools, Workshop Welding Wire Tool Kits, Workshop Welding Wire Manufacturers, Workshop Welding Wire Suppliers from India, China, Kenya



**Laboratory Instrument India**

