




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
FWD Diesel Engine With A Gearbox Cutaway

**Product Code :**  
LBNY-0005-121000450



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Front Wheel Drive (FWD) Diesel Engine With A Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

**Technical Specification :**

4 stroke engine  
4 in-line cylinders  
Distribution through a toothed belt  
Thermostatic valve  
Alternator  
Displacement: 1300 cc  
Gearbox: 5 forward speeds + reverse and differential  
Single-plate clutch with diaphragm  
Rotary injection pump  
Indirect injection  
Overhead camshaft (OHC)  
Power: 45 hp.

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 130 x 86 x 100(h) cm

---

Net Weight: 140 kg  
Gross Weight: 200 kg  
Equipment Layout: On Stand with Wheels



**Laboratory Instrument India**