




Email : sales@laboratoryinstrumentindia.com

Product Name :
FWD Diesel Engine With A Gearbox Cutaway

Product Code :
LBNY-0005-121000450



Description :

Many parts have been chromium-plated and galvanized for a longer life. Front Wheel Drive (FWD) Diesel Engine With A Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

4 stroke engine
4 in-line cylinders
Distribution through a toothed belt
Thermostatic valve
Alternator
Displacement: 1300 cc
Gearbox: 5 forward speeds + reverse and differential
Single-plate clutch with diaphragm
Rotary injection pump
Indirect injection
Overhead camshaft (OHC)
Power: 45 hp.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 130 x 86 x 100(h) cm

Net Weight: 140 kg
Gross Weight: 200 kg
Equipment Layout: On Stand with Wheels



Laboratory Instrument India