




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
HGV Diesel Engine with Electronically Controlled Pump  
Injectors Cutaway

**Product Code :**  
LBNY-0005-121000454



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. HGV (Heavy Goods Vehicle )Diesel Engine with Electronically Controlled Pump Injectors Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

**Technical Specification :**

4 stroke  
4 valves per cylinder  
6 in-line cylinders  
Turbo-compressor  
Pump injectors electronically controlled  
Displacement: 7790/10380 cc  
Maximum power 310/450hp  
Water cooling  
Pre-heating device.

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 200 x 105 x 150(h) cm  
Net Weight: 950 kg  
Gross Weight: 1200 kg

---

Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**