




Email : sales@laboratoryinstrumentindia.com

Product Name :
HGV Diesel Engine with Electronically Controlled Pump
Injectors Cutaway

Product Code :
LBNY-0005-121000454



Description :

Many parts have been chromium-plated and galvanized for a longer life. HGV (Heavy Goods Vehicle)Diesel Engine with Electronically Controlled Pump Injectors Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

4 stroke
4 valves per cylinder
6 in-line cylinders
Turbo-compressor
Pump injectors electronically controlled
Displacement: 7790/10380 cc
Maximum power 310/450hp
Water cooling
Pre-heating device.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 200 x 105 x 150(h) cm
Net Weight: 950 kg
Gross Weight: 1200 kg

Equipment Layout: On Stand with Wheels.



Laboratory Instrument India