




Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbo Diesel Engine For Car And Lorry With A Gearbox
Cutaway

Product Code :
LBNY-0005-121000456



Description :

Many parts have been chromium-plated and galvanized for a longer life. Turbo Diesel Engine For Car And Lorry With A Gearbox Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc.

Technical Specification :

4 in-line cylinders
4 stroke engine
Overhead camshaft (OHC)
Distribution through a toothed belt;
Gearbox: 5 forward speeds + reverse
Water cooling
Alternator
Thermostatic valve
Displacement: 2000/2500 cc
Indirect injection
Rotary injection pump

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 85 x 160 x 100(h) cm
Net Weight: 220 kg

Gross Weight: 280 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India