



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Spares

Product Code :
LBNY-0004-740014



Description :

Workshop Spares

Technical Specification :

Spares

MAGS spare welding tips, shrouds, torches and cutting nozzles.

TAGS Tungstens, collets, defusers, back caps ceramic shrouds.

Spare MAGS torch and TAGS torch along with MMA torches and return cables and clamps.

Spare lenses for welding helmets clear and tinted

Spare burning and grinding visors

Workshop Spares, Workshop Spares Equipments, Workshop Spares Tools, Workshop Spares Tool Kits,
Workshop Spares Manufacturers, Workshop Spares Suppliers from India, China, Kenya



Laboratory Instrument India