



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tracked Diesel Tractor FIAT Cutaway

Product Code :
LBNY-0005-121000460



Description :

Many parts have been chromium-plated and galvanized for a longer life. Tracked Diesel Tractor FIAT Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, lubricating circuits, cross-sections, cooling system etc.

Technical Specification :

Indirect injection
Water cooling
Overhead valves
4 stroke 4 cylinder engine
Displacement: 2000 cc
Gearbox: 4 forward speeds + reverse
In-line injection pump
Multiple-plate steering clutch.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 230 x 132 x 180(h) cm
Net Weight: 910 kg
Gross Weight: 1200 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India