




Email : sales@laboratoryinstrumentindia.com

Product Name :
HGV Turbo Diesel Engine With Intercooler And A Gearbox
Cutaway

Product Code :
LBNY-0005-121000461



Description :

Many parts have been chromium-plated and galvanized for a longer life. HGV (Heavy Goods Vehicle) Turbo Diesel Engine With Intercooler And A Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, lubricating circuits, cross-sections, cooling system etc.

Technical Specification :

4 stroke engine
6 in-line cylinders
Intercooler air-air
Fuel fed by turbo-supercharger and intercooler
Gearbox: 4 speeds and reverse
Camshaft in the crankcase
Spring single-plate clutch
Geared distribution
Displacement: 5200 cc
Direct injection
In-line injection pump with mechanical governor

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 90 x 180 x 170(h) cm

Net Weight: 450 kg
Gross Weight: 560 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India