



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle Cutaway

Product Code :
LBNY-0005-121000463



Description :

Many parts have been chromium-plated and galvanized for a longer life. Motorcycle Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, fuel system, cross-sections, cooling system etc.

Technical Specification :

Operated: Electrically (Motor 230 V)
Dimensions: 200 x 70 x 130(h) cm
Net Weight: 160 kg
Gross Weight: 230 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India