



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle Cutaway

**Product Code :**  
LBNY-0005-121000463



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Motorcycle Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, fuel system, cross-sections, cooling system etc.

**Technical Specification :**

Operated: Electrically (Motor 230 V)  
Dimensions: 200 x 70 x 130(h) cm  
Net Weight: 160 kg  
Gross Weight: 230 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**