



Email : sales@laboratoryinstrumentindia.com

Product Name :
Marine Inboard Diesel Engine With Inverter Cutaway

Product Code :
LBNY-0005-121000464



Description :

Many parts have been chromium-plated and galvanized for a longer life. Marine Inboard Diesel Engine With Inverter Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, fuel system, cross-sections, cooling system etc.

Technical Specification :

Operated: Electrically (Motor 230 V)
4 or 6-cylinder Diesel Engine
Dimensions and Weight depend on engine type (4 or 6 cylinders).
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India