




Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle Two Stroke Petrol Engine Model Cutaway

Product Code :
LBNY-0005-121000468



Description :

Many parts have been chromium-plated and galvanized for a longer life. Motorcycle Two-stroke Petrol Engine Model Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

Plug point and magnet flywheel ignition
Piston displacement 48 cc
Air cooling
Box carburettor

Other Specifications:

Operated: Manually (Through a Crank Handle)
Dimensions: 24 x 42 x 28(h) cm
Net Weight: 6 kg
Gross Weight: 10 kg
Equipment Layout: On Base.

Laboratory Instrument India

