




Email : sales@laboratoryinstrumentindia.com

Product Name :
Marine Inboard Diesel Engine Without Inverter Cutaway

Product Code :
LBNY-0005-121000472



Description :

Many parts have been chromium-plated and galvanized for a longer life. Marine Inboard Diesel Engine Without Inverter Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

Operated: Electrically (Motor 230 V)
Dimensions and Weight depend on engine type (6 or 8 cylinders)
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India