



Email : sales@laboratoryinstrumentindia.com

Product Name :
FWD Chassis With Turbo Diesel Engine Fully Functional
Cutaway

Product Code :
LBNY-0005-121000473



Description :

All the mechanical parts are carefully revised and assembled on a chassis. FWD Chassis with Turbo Diesel Engine Fully Functional Cutaway trainer is realized from a medium displacement vehicle which is particularly useful for the study of the main components, the repair operations and the diagnostic, by means of the OBD socket. This chassis allows an easy display and the possibility of assemble and disassemble each part.

Technical Specification :

- 4 cylinder 4 stroke Turbo Diesel Engine
- Electronic ignition
- Displacement: 1700 cc
- Gearbox: 5 forward gears and reverse
- VE injection pump
- Depression servo-brake
- Front wheel drive
- Silencer
- Rear drum brakes
- McPherson front suspension
- Front disc brakes
- Independent wheels rear suspension with oscillating arms
- Battery
- Radiator with electric fan
- Hand brake

Tank
Dashboard

Other Specifications:

Operated: Fully Functional

Dimensions: 300 x 150 x 140(h) cm

Net Weight: 450 kg

Gross Weight: 730 kg

Equipment Layout: On Stand with Wheels.



Laboratory Instrument India