



Email : sales@laboratoryinstrumentindia.com

Product Name :

FWD Chassis With Turbo Diesel Engine Cutaway

Product Code :

LBNY-0005-121000474



Description :

Many parts have been chromium-plated and galvanized for a longer life. FWD Chassis With Turbo Diesel Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

Technical Specification :

Fiat chassis with front drive (engine transversally mounted)
Diesel engine displacement: 1700 cc
Gearbox: 5 forward speed+reverse+differential
Hydraulic power steering with double-jointed steering column
Radiator with electric fan
4 cylinders
Indirect injection with pre-chamber
VE rotary injection pump
Turbo-compressor with waste-gate valve
Independent wheels McPherson front suspension with oscillating arms
Rear independent suspension with coil spring and gas shock absorbers
Front disc brake
Rear drum brake

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 220 x 160 x 115(h) cm

Net Weight: 290 kg
Gross Weight: 400 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India