



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Outboard Marine Reverser Cutaway

**Product Code :**  
LBNY-0005-121000475



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Outboard Marine Reverser Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, fuel system, cross-sections, cooling system etc.

**Technical Specification :**

Operated: Manually (Through a Crank Handle)  
Dimensions: 50 x 40 x 75(h) cm  
Net Weight: 8 kg  
Gross Weight: 15 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**