



Email : sales@laboratoryinstrumentindia.com

Product Name :
Outboard Marine Reverser Cutaway

Product Code :
LBNY-0005-121000475



Description :

Many parts have been chromium-plated and galvanized for a longer life. Outboard Marine Reverser Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, fuel system, cross-sections, cooling system etc.

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 50 x 40 x 75(h) cm
Net Weight: 8 kg
Gross Weight: 15 kg
Equipment Layout: On Base.



Laboratory Instrument India