




Email : sales@laboratoryinstrumentindia.com

Product Name :
Common Rail Turbo Diesel Engine Cutaway

Product Code :
LBNY-0005-121000481



Description :

Many parts have been chromium-plated and galvanized for a longer life. Common Rail Turbo Diesel Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

4 stroke engine
4 in-line cylinders
Power: 115 hp. At 4000 RPM
2 valves per cylinder
Displacement: 1900 cc
Alternator-oil filter-oil pump
Common rail-type direct injection with electro- injectors;
Overhead camshaft (OHC) with timing belt
Turbo-supercharger.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 90 x 100 x 120(h) cm
Net Weight: 150 kg
Gross Weight: 210 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India