




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Common Rail Turbo Diesel Engine Cutaway

**Product Code :**  
LBNY-0005-121000481



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Common Rail Turbo Diesel Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

**Technical Specification :**

4 stroke engine  
4 in-line cylinders  
Power: 115 hp. At 4000 RPM  
2 valves per cylinder  
Displacement: 1900 cc  
Alternator-oil filter-oil pump  
Common rail-type direct injection with electro- injectors;  
Overhead camshaft (OHC) with timing belt  
Turbo-supercharger.

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 90 x 100 x 120(h) cm  
Net Weight: 150 kg  
Gross Weight: 210 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**