




Email : sales@laboratoryinstrumentindia.com

Product Name :
Direct Fuel Injection Turbo Diesel Engine For Car And Lorry
Cutaway

Product Code :
LBNY-0005-121000482



Description :

Direct Fuel Injection Turbo Diesel Engine For Car And Lorry Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc. Many parts have been chromium-plated and galvanized for a longer life.

Technical Specification :

4 stroke engine
4 in-line cylinders
Feeding by turbo-supercharger
Rotary injection pump
Direct/Indirect injection
Overhead camshaft (OHC)
Displacement: 2500 cc
Distribution through a toothed belt
Alternator
Water cooling
Thermostatic valve
Intercooler water-oil.

Other Specifications:

Operated: Electrically (Motor 230 V)

Dimensions: 90 x 100 x 120(h) cm

Net Weight: 200 kg

Gross Weight: 260 kg

Equipment Layout: On Stand with Wheels.



Laboratory Instrument India