




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Direct Fuel Injection Turbo Diesel Engine For Car And Lorry  
Cutaway

**Product Code :**  
LBNY-0005-121000482



**Description :**

Direct Fuel Injection Turbo Diesel Engine For Car And Lorry Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc. Many parts have been chromium-plated and galvanized for a longer life.

**Technical Specification :**

4 stroke engine  
4 in-line cylinders  
Feeding by turbo-supercharger  
Rotary injection pump  
Direct/Indirect injection  
Overhead camshaft (OHC)  
Displacement: 2500 cc  
Distribution through a toothed belt  
Alternator  
Water cooling  
Thermostatic valve  
Intercooler water-oil.

**Other Specifications:**

Operated: Electrically (Motor 230 V)

---

Dimensions: 90 x 100 x 120(h) cm

Net Weight: 200 kg

Gross Weight: 260 kg

Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**