




Email : sales@laboratoryinstrumentindia.com

Product Name :

RWD Turbo Diesel Engine with a Gearbox Cutaway

Product Code :

LBNY-0005-121000483



Description :

Many parts have been chromium-plated and galvanized for a longer life. Rear Wheel Drive (RWD) Turbo Diesel Engine with a Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

Technical Specification :

Displacement: 1300/1700 cc
Indirect injection
4 stroke engine
4 in-line cylinders
Gearbox: 4 forward speeds + reverse
Single-plate clutch with diaphragm
Water cooling
Feeding by turbo supercharger;
Rotary injection pump
Alternator
Thermostatic valve
Overhead camshaft (OHC)
Distribution through a toothed belt.

Other Specifications:

Operated: Electrically (Motor 230 V)

Dimensions: 60 x 160 x 100(h) cm

Net Weight: 150 kg

Gross Weight: 210 kg

Equipment Layout: On Stand with Wheels.



Laboratory Instrument India