




Email : sales@laboratoryinstrumentindia.com

Product Name :
RWD Chassis With DOHC Multipoint EFI Petrol Engine
Cutaway

Product Code :
LBNY-0005-121000486



Description :

Many parts have been chromium-plated and galvanized for a longer life. RWD Chassis with DOHC Multipoint EFI Petrol Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

4 stroke 4 cylinder petrol engine
Gearbox: 5 forward speeds + reverse,
Displacement: 2000 cc
Differential with hypoid crown wheel and pinion;
Twin overhead camshaft driven by a toothed belt;
Front disc brakes and rear drum brakes
Rack steering box
Electronic ignition
Dual braking circuit
Rear leaf spring suspension.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 145 x 220 x 100(h) cm
Net Weight: 300 kg

Gross Weight: 420 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India