




Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Cylinder Four Stroke Diesel Engine Air Cooled
Cutaway

Product Code :
LBNY-0005-121000487



Description :

Many parts have been chromium-plated and galvanized for a longer life. Single Cylinder Four Stroke Diesel Engine Air Cooled Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

Power 4 hp Direct Injection
Camshaft in the Crankcase
Displacement 210 cc
Overhead Camshafts
RPM Regulator
Injector
Silencer
Toroidal Oil Pump
Injection Pump

Other Specifications:

Operated: Manually
Dimensions: 40 x 40 x 50(h) cm
Net Weight: 25 kg
Gross Weight: 35 kg

Equipment Layout: On Base.



Laboratory Instrument India