




Email : sales@laboratoryinstrumentindia.com

Product Name :
Carburettor Petrol Engine with Gearbox Cutaway

Product Code :
LBNY-0005-121000490



Description :

Many parts have been chromium-plated and galvanized for a longer life. Carburettor Petrol Engine (FIAT) with Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc.

Technical Specification :

4 Cylinders
Coil ignition
Mechanical petrol pump
Gearbox 4 speeds and reverse
Displacement: 1200/2000 cc
In-line overhead valves
Dry single-plate clutch
Camshaft in the crankcase, belt/chain drive.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 60 x 155 x 90(h) cm
Net Weight: 140 kg
Gross Weight: 185 kg.
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India