




Email : sales@laboratoryinstrumentindia.com

Product Name :
Twin rotor Wankel Petrol Engine Cutaway

Product Code :
LBNY-0005-121000493



Description :

Many parts have been chromium-plated and galvanized for a longer life. Twin-rotor Wankel Petrol Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

Operated: Manually
Dimensions: 70 x 70 x 100(h) cm
Net Weight: 95 kg
Gross Weight: 140 kg.
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India