




Email : sales@laboratoryinstrumentindia.com

Product Name :
Carburettor Twin Shaft Petrol Engine With Gearbox
Cutaway

Product Code :
LBNY-0005-121000495



Description :

Many parts have been chromium-plated and galvanized for a longer life. Carburettor Twin-Shaft Petrol Engine With Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc.

Technical Specification :

4 Cylinders
Displacement: 1600/2000 cc
Gearbox: 5 speeds and reverse
2 overhead camshaft driven by a toothed belt
Twin-carburettor
Overhead valves with V-arrangement
Coil ignition
Alternator
Dry single-plate clutch

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 76 x 140 x 100(h) cm
Net Weight: 150 kg
Gross Weight: 220 kg

Equipment Layout: On Stand with Wheels.



Laboratory Instrument India