




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
RWD L Jetronic EFI Petrol Engine with Gearbox Cutaway

**Product Code :**  
LBNY-0005-121000500



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. RWD L Jetronic EFI Petrol Engine with Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc.

**Technical Specification :**

4 in-line cylinders  
Gearbox: 5 speeds and reverse  
Max power: 140 hp.  
Displacement: 2000 cc  
2 overhead camshafts  
L-jetronic electronic injection  
Electronic ignition

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 80 x 140 x 110(h) cm  
Net Weight: 140 kg  
Gross Weight: 200 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**