




Email : sales@laboratoryinstrumentindia.com

Product Name :
RWD L Jetronic EFI Petrol Engine with Gearbox Cutaway

Product Code :
LBNY-0005-121000500



Description :

Many parts have been chromium-plated and galvanized for a longer life. RWD L Jetronic EFI Petrol Engine with Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, fuel system, lubricating circuits, cooling system etc.

Technical Specification :

4 in-line cylinders
Gearbox: 5 speeds and reverse
Max power: 140 hp.
Displacement: 2000 cc
2 overhead camshafts
L-jetronic electronic injection
Electronic ignition

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 80 x 140 x 110(h) cm
Net Weight: 140 kg
Gross Weight: 200 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India