



Email : sales@laboratoryinstrumentindia.com

Product Name :
Petrol Engine Cutaway

Product Code :
LBNY-0005-121000503



Description :

Many parts have been chromium-plated and galvanized for a longer life. Petrol Engine Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

Technical Specification :

4 cylinders
Displacement: 1000-1300 cc
Multi-point electronic injection with throttle
12V alternator
Thermostatic valve
DOHC - double overhead camshaft
Roller chain
4 valves per cylinder.

Other Specifications:

Operated: Manually (Through a Crank Handle)
Dimensions: 70 x 90 x 100(h) cm
Net Weight: 60 kg
Gross Weight: 110 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India