




Email : sales@laboratoryinstrumentindia.com

Product Name :
Multipoint Electronic Fuel Injection DOHC Petrol Engine
Gearbox Cutaway

Product Code :
LBNY-0005-121000505



Description :

Many parts have been chromium-plated and galvanized for a longer life. Multipoint Electronic Fuel Injection DOHC Petrol Engine Gearbox Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, lubricating circuits, cross-sections, fuel system, cooling system etc.

Technical Specification :

4 valves per cylinder
6 in-line cylinders
Gearbox: 5 forward speeds and reverse
Displacement: 2000/2500 cc
Double overhead camshaft (DOHC)
Distribution by means of a roller chain

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 140 x 80 x 100(h) cm
Net Weight: 180 kg
Gross Weight: 265 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India