




Email : sales@laboratoryinstrumentindia.com

Product Name :
Petrol Engine Ferrari Cutaway

Product Code :
LBNY-0005-121000506



Description :

Many parts have been chromium-plated and galvanized for a longer life. Petrol Engine (Ferrari) Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

Technical Specification :

4 variable timing devices on the camshaft
Multi-point electronic injection
Geared distribution
12 V cylinders
4 valves per cylinder
Double disc clutch
Water pump
DOHC overhead camshaft
3 oil pumps.

Other Specifications:

Operated: Manually (Through a Crank Handle)
Dimensions: 80 x 100 x 130(h) cm
Net Weight: 180 kg
Gross Weight: 250 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India