




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Petrol Engine Ferrari Cutaway

**Product Code :**  
LBNY-0005-121000506



#### **Description :**

Many parts have been chromium-plated and galvanized for a longer life. Petrol Engine (Ferrari) Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

#### **Technical Specification :**

4 variable timing devices on the camshaft  
Multi-point electronic injection  
Geared distribution  
12 V cylinders  
4 valves per cylinder  
Double disc clutch  
Water pump  
DOHC overhead camshaft  
3 oil pumps.

#### **Other Specifications:**

Operated: Manually (Through a Crank Handle)  
Dimensions: 80 x 100 x 130(h) cm  
Net Weight: 180 kg  
Gross Weight: 250 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**